

GDAP1L1 Antibody



PACO09440

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

GDAP1L1

Uniprot

Q96MZ0

Synonyms:

ganglioside-induced differentiation-associated protein 1-like 1; GDAP1L1;
DKFZp761K228; dJ881L22.1; dJ995J12.1.1 ;

Immunogen:

Human GDAP1L1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A