

# GDE1 Antibody



PACO09441

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

GDE1

**Uniprot**

Q9NZC3

**Synonyms:**

glycerophosphodiester phosphodiesterase 1; GDE1;363E6.2; MIR16 ;

**Immunogen:**

Human GDE1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A