

GGT5 Antibody



PACO09460

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

GGT5

Uniprot

P36269

Synonyms:

gamma-glutamyltransferase 5; GGT5; DKFZp566O011; FLJ92733; GGT-REL; GGTLA1 ;

Immunogen:

Human GGT5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A