

# GH2 Antibody



PACO09462

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

GH2

**Uniprot**

P01242

**Synonyms:**

growth hormone 2; GH2; GH-V; GHL; GHV; hGH-V ;

**Immunogen:**

Human GH2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A