

PACO09463

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

GHITM

**Uniprot**

Q9H3K2

**Synonyms:**

growth hormone inducible transmembrane protein; GHITM; DERP2; DKFZp566C0746; FLJ26584; HSPC282; MICS1; My021; PTD010; TMBIM5 ;

**Immunogen:**

Human GHITM.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A