

PACO09464

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

GHRL

Uniprot

Q9UBU3

Synonyms:

ghrelin/obestatin prepropeptide; GHRL; MTLRP; obestatin ;

Immunogen:

Human GHRL.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A