

GIMAP4 Antibody



PACO09466

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GIMAP4

Uniprot

Q9NUV9

Synonyms:

GTPase, IMAP family member 4; GIMAP4; FLJ11110; HIMAP4; IAN1; IMAP4; MSTP062; hIAN1 ;

Immunogen:

Human GIMAP4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A