

GINS2 Antibody



PACO09469

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

GINS2

Uniprot

Q9Y248

Synonyms:

GINS complex subunit 2 (Psf2 homolog); GINS2; HSPC037; PSF2; Pfs2 ;

Immunogen:

Human GINS2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A