

PACO09470

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

GINS3

Uniprot

Q9BRX5

Synonyms:

GINS complex subunit 3 (Psf3 homolog); GINS3; FLJ13912; PSF3 ;

Immunogen:

Human GINS3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A