

GK5 Antibody



PACO09476

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

GK5

Uniprot

Q6ZS86

Synonyms:

glycerol kinase 5 (putative); GK5; FLJ33582; FLJ45739; MGC40579 ;

Immunogen:

Human GK5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A