

GLI4 Antibody



PACO09480

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GLI4

Uniprot

P10075

Synonyms:

GLI family zinc finger 4; GLI4; HKR4 ;

Immunogen:

Human GLI4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A