

GLIPR1L1 Antibody



PACO09482

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

GLIPR1L1

Uniprot

Q6UWM5

Synonyms:

GLI pathogenesis-related 1 like 1; GLIPR1L1; ALKN2972; MGC26856; PRO7434 ;

Immunogen:

Human GLIPR1L1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A