

GLO1 Antibody



PACO09483

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

GLO1

Uniprot

Q04760

Synonyms:

glyoxalase I; GLO1; GLOD1; GLYI ;

Immunogen:

Human GLO1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A