

GLRA3 Antibody



PACO09486

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

GLRA3

Uniprot

O75311

Synonyms:

glycine receptor, alpha 3;

Immunogen:

Human GLRA3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; Avoid freeze / thaw cycles.

Product Images

N/A

N/A