

GM2A Antibody



PACO09497

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

GM2A

Uniprot

P17900

Synonyms:

GM2 ganglioside activator; GM2A; SAP-3 ;

Immunogen:

Human GM2A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A