

GNAO1 Antibody



PACO09510

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

GNAO1

Uniprot

P09471

Synonyms:

guanine nucleotide binding protein (G protein), alpha activating activity polypeptide O; GNAO1; DKFZp686O0962; G-ALPHA-o; GNAO ;

Immunogen:

Human GNAO1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A