

GNB5 Antibody



PACO09518

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

GNB5

Uniprot

O14775

Synonyms:

guanine nucleotide binding protein (G protein), beta 5; GNB5; FLJ37457; FLJ43714; GB5 ;

Immunogen:

Human GNB5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A