

GOT2 Antibody



PACO09540

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

GOT2

Uniprot

P00505

Synonyms:

glutamic-oxaloacetic transaminase 2, mitochondrial (aspartate aminotransferase 2);
GOT2; FLJ40994; KAT4; KATIV; mitAAT ;

Immunogen:

Human GOT2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A