

GPBP1L1 Antibody



PACO09544

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GPBP1L1

Uniprot

Q9HC44

Synonyms:

GC-rich promoter binding protein 1-like 1; GPBP1L1; RP11-767N6.1; SP192 ;

Immunogen:

Human GPBP1L1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A