

GPR137 Antibody



PACO09559

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

GPR137

Uniprot

Q96N19

Synonyms:

G protein-coupled receptor 137; GPR137; C11orf4; GPR137A; TM7SF1L1 ;

Immunogen:

Human GPR137.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A