

GPRASP2 Antibody



PACO09570

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

GPRASP2

Uniprot

Q96D09

Synonyms:

G protein-coupled receptor associated sorting protein 2; GPRASP2; FLJ35662; FLJ37327; GASP2 ;

Immunogen:

Human GPRASP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A