

GPSM2 Antibody



PACO09574

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

GPSM2

Uniprot

P81274

Synonyms:

G-protein signaling modulator 2 (AGS3-like, C. elegans); GPSM2; LGN; Pins ;

Immunogen:

Human GPSM2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A