

# GRB2 Antibody



PACO09584

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

GRB2

**Uniprot**

P62993

**Synonyms:**

growth factor receptor-bound protein 2; GRB2; ASH; EGFRBP-GRB2; Grb3-3; MST084; MSTP084 ;

**Immunogen:**

Human GRB2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A