

# GRID1 Antibody



PACO09588

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

GRID1

**Uniprot**

Q9ULK0

**Synonyms:**

glutamate receptor, ionotropic, delta 1; GRID1; KIAA1220 ;

**Immunogen:**

Human GRID1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A