

GRK6 Antibody



PACO09595

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

GRK6

Uniprot

P43250

Synonyms:

G protein-coupled receptor kinase 6; GRK6; FLJ32135; GPRK6 ;

Immunogen:

Human GRK6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A