

GSTM3 Antibody



PACO09616

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

GSTM3

Uniprot

P21266

Synonyms:

glutathione S-transferase mu 3 (brain); GSTM3; GST5; GSTB; GSTM3-3; GTM3; MGC3310; MGC3704 ;

Immunogen:

Human GSTM3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A