

GSX1 Antibody



PACO09625

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

GSX1

Uniprot

Q9H4S2

Synonyms:

GS homeobox 1; GSX1; GSH1; Gsh-1 ;

Immunogen:

Human GSX1.

Storage:

0.1M NaHCO₃, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A