

GTF3C5 Antibody



PACO09637

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GTF3C5

Uniprot

Q9Y5Q8

Synonyms:

general transcription factor IIIC, polypeptide 5, 63kDa; GTF3C5; FLJ20857; TFIIIC63; TFIIICepsilon; TFIIIC2-63 ;

Immunogen:

Human GTF3C5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A