

GTPBP3 Antibody



PACO09640

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

GTPBP3

Uniprot

Q969Y2

Synonyms:

GTP binding protein 3 (mitochondrial); GTPBP3; FLJ14700; GTPBG3; MSS1; MTGP1; THDF1 ;

Immunogen:

Human GTPBP3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A