

MTG2 Antibody



PACO09642

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MTG2

Uniprot

Q9H4K7

Synonyms:

GTP binding protein 5 (putative); GTPBP5; FLJ10741; MGC29512; ObgH1; dJ1005F21.2 ;

Immunogen:

Human GTPBP5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A