

GULP1 Antibody



PACO09649

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GULP1

Uniprot

Q9UBP9

Synonyms:

GULP, engulfment adaptor PTB domain containing 1; GULP1; CED-6; CED6; FLJ31156; GULP ;

Immunogen:

Human GULP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A