
Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

GYPA

Uniprot

P02724

Synonyms:

glycophorin A (MNS blood group); GYPA; CD235a; GPA; GP_Erik; GPSAT; Gp^{MI}III; HGp^{MI}III; HGp^{MI}IV; HGp^{MI}X; HGp^{MI}XI; HGpSta(C); MN; MNS ;

Immunogen:

Human GYPA.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; Avoid freeze / thaw cycles.

Product Images

N/A

N/A