

GZMM Antibody



PACO09656

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GZMM

Uniprot

P51124

Synonyms:

granzyme M (lymphocyte met-ase 1); GZMM; LMET1; MET1 ;

Immunogen:

Human GZMM.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A