

H6PD Antibody



PACO09663

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

H6PD

Uniprot

O95479

Synonyms:

hexose-6-phosphate dehydrogenase (glucose 1-dehydrogenase); H6PD; DKFZp686A01246; G6PDH; GDH; MGC87643 ;

Immunogen:

Human H6PD.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A