

HADH Antibody



PACO09667

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

HADH

Uniprot

Q16836

Synonyms:

hydroxyacyl-Coenzyme A dehydrogenase; HADH; HAD; HADH1; HADHSC; HHF4; M/SCHAD; MGC8392; SCHAD ;

Immunogen:

Human HADH.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A