

HARS Antibody



PACO09671

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

HARS

Uniprot

P12081

Synonyms:

histidyl-tRNA synthetase; HARS; FLJ20491; HRS ;

Immunogen:

Human HARS.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A