

HBB Antibody



PACO09677

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

HBB

Uniprot

P68871

Synonyms:

hemoglobin, beta; HBB; CD113t-C; beta-globin ;

Immunogen:

Human HBB.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A