

HCFC2 Antibody



PACO09686

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

HCFC2

Uniprot

Q9Y5Z7

Synonyms:

host cell factor C2; HCFC2; FLJ94012; HCF-2; HCF2 ;

Immunogen:

Human HCFC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A