

HCN3 Antibody



PACO09689

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

HCN3

Uniprot

Q9P1Z3

Synonyms:

hyperpolarization activated cyclic nucleotide-gated potassium channel 3; HCN3; KIAA1535; MGC131493 ;

Immunogen:

Human HCN3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A