

# HDAC1 Antibody



PACO09691

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

HDAC1

**Uniprot**

Q13547

**Synonyms:**

histone deacetylase 1; HDAC1; DKFZp686H12203; GON-10; HD1; RPD3; RPD3L1 ;

**Immunogen:**

Human HDAC1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A