

HDGF Antibody



PACO09698

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

HDGF

Uniprot

P51858

Synonyms:

hepatoma-derived growth factor (high-mobility group protein 1-like); HDGF; DKFZp686J1764; FLJ96580; HMG1L2 ;

Immunogen:

Human HDGF.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A