

HDLBP Antibody



PACO09700

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

HDLBP

Uniprot

Q00341

Synonyms:

high density lipoprotein binding protein; HDLBP; FLJ16432; HBP; PRO2900; VGL ;

Immunogen:

Human HDLBP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A