

HEBP2 Antibody



PACO09702

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

HEBP2

Uniprot

Q9Y5Z4

Synonyms:

heme binding protein 2; HEBP2; C6ORF34B; C6orf34; KIAA1244; PP23; SOUL ;

Immunogen:

Human HEBP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A