

HECTD1 Antibody



PACO09703

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

HECTD1

Uniprot

Q9ULT8

Synonyms:

HECT domain containing 1; HECTD1; FLJ38315; KIAA1131 ; ; EULIR; HECD1

Immunogen:

Human HECTD1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A