

# HERPUD1 Antibody



PACO09710

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

HERPUD1

**Uniprot**

Q15011

**Synonyms:**

homocysteine-inducible, endoplasmic reticulum stress-inducible, ubiquitin-like domain member 1; HERPUD1; HERP; KIAA0025; Mif1; SUP ;

**Immunogen:**

Human HERPUD1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A