

HEXIM1 Antibody



PACO09713

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

HEXIM1

Uniprot

O94992

Synonyms:

hexamethylene bis-acetamide inducible 1; HEXIM1; CLP1; EDG1; FLJ13562; HIS1; MAQ1 ;

Immunogen:

Human HEXIM1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A