

HFE Antibody



PACO09716

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

HFE

Uniprot

Q30201

Synonyms:

hemochromatosis; HFE; HFE1; HH; HLA-H; MGC103790; MVCD7; dJ221C16.10.1 ;

Immunogen:

Human HFE.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A