

HIF1A Antibody



PACO09723

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HIF1A

Uniprot

Q16665

Synonyms:

hypoxia inducible factor 1, alpha subunit (basic helix-loop-helix transcription factor); HIF1A; HIF-1alpha; HIF1; HIF1-ALPHA; MOP1; PASD8; bHLHe78 ;

Immunogen:

Human HIF1A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A