

# HIST3H2A Antibody



PACO09732

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

HIST3H2A

**Uniprot**

Q7L7L0

**Synonyms:**

histone cluster 3, H2a; HIST3H2A; MGC3165 ;

**Immunogen:**

Human HIST3H2A.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A