

HM13 Antibody



PACO09739

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HM13

Uniprot

Q8TCT9

Synonyms:

histocompatibility (minor) 13; HM13; H13; IMP1; IMPAS; MSTP086; PSENL3; PSL3; SPP; dJ324O17.1 ;

Immunogen:

Human HM13.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A